### REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N 302173

RECALIBRATION DATE 3-30-0/

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34-00	
37.00	37.00	37.00	

Tested by D. H. H.
Date 3-3-00
Approved
Date

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N 09/802

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30,00	<b>3</b> 0.60	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by A. W.
Date <u>3-3-00</u>
Approved
Date

## REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>09/800</u>

RECALIBRATION DATE 3-3/-0/

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	32.00	

Tested by A. W.
Date 3-30-00
Approved
Date

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>09/793</u>

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by A. H.
Date 3-3-00
Approved
Date

## REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>09/794</u>

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	3700	

Tested by A. W.
Date 3-3-00
Approved
Date

## REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>09/795</u>

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by A. V.
Date 3-3-00
Approved
Date

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>091796</u>

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30,00	
34.00	34:00	34.00	
37.00	37.00	37.00	

Tested by A.W.		
Date	3-3-00	
Approved		
Date		

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>091797</u>

RECALIBRATION DATE 3-30-0/

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by A. W.
Date 3-3-60
Approved
Date

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>09/798</u>

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by N.W.	_
Date 3-3-00	-
Approved	_
Date	

REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER  $29.5^{\circ}$  C  $-38.5^{\circ}$  C

GUTH P/N GU-4300-DT

S/N 09/799

RECALIBRATION DATE 3-3/-0/

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by 💋	D.A. W.
Date_	3-3-00
Approved_	
Date	

REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER  $29.5^{\circ}$  C  $-38.5^{\circ}$  C

GUTH P/N GU-4300-DT

S/N 09/80/

RECALIBRATION DATE 3-30-01

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	
34.00	34.00	34.00	
37.00	37.00	37.00	

Tested by A. A.
Date 3-3-00
Approved
Date

# REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER $29.5^{\circ}$ C $-38.5^{\circ}$ C

GUTH P/N GU-4300-DT

S/N <u>082708</u>

RECALIBRATION DATE <u>1/18/03</u>

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 300905	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	1/18/02
34.00	34.00	34.00	1/18/02
37.00	37.00	37.00	1/18/02

Tested by _	131
Date_	1/18/2002
Approved_	
Date	

### REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER 29.5° C – 38.5° C

GUTH P/N GU-4300-DT

S/N <u>082708</u>

RECALIBRATION DATE \_\_/-//-2002\_

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 300905	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00	30.00	1-11-01
34.00	34.00	34.00	1-11-01
37.00	37 60	37.00	1-11-01

Tested by	DAH.
Date_	1-11-2001
Approved_	
Date	

REPORT OF CALIBRATION GUTH-EUTECHNICS DIGITAL THERMOMETER  $29.5^{\circ}$  C  $-38.5^{\circ}$  C

GUTH P/N GU-4300-DT

S/N <u>082708</u>

RECALIBRATION DATE 12/2/99

The Digital Thermometer calibration was verified at 30.00° C, 34.00° C, and 37.00° C by comparison to temperature standards traceable to the National Institute of Standards and Technology.

BATH TEMPERATURE (°C)	MASTER DIGITAL THERMOMETER S/N 074857	DIGITAL THERMOMETER READING	TEST DATE
30.00	30.00°C	30.00°C	12/99
34.00	34.00°C	34.00%	14/2/99
37.00	37.00°C	37.00°C	12/2/99

Tested by 🚄

Date 13/3/

Approved

Date 12/2/99